

# HaloGX Pro™

In response to the coronavirus, more and more spas are creating smaller, more personalized salt room and cave experiences. We have also adapted and updated our technology to introduce the HaloGX Pro. The HaloGX Pro produces the smallest particles in the industry and uses the newest hybrid methods, incorporating both mill technology and cyclone technology.

## FEATURES

- Great for rooms up to 150+ square feet
- Engineered with new compressed air technology for the smallest, micron-sized salt particles in the industry
- Uses less salt per session than standard Halogenerators but with the same efficacy
- Very quiet machine
- Simple installation (2 screws and a plug)
- Very easy to operate and clean
- Includes "turbo" feature to increase salt levels if desired
- Low power consumption 50W, 110/230V, 50/60Hz
- LIFETIME Warranty
- Perfect size: 26.1 (H) x 8.3 (W) x 6.6 (D) (inches)
- FREE Consultation on build-out

## LIST PRICE

- **\$4,900**

*\*For all customers, shipping cost is not included. For Non-US customers, you are responsible for customs/duties when applicable*



# HaloFX™

Our unique HaloFX™ Halogenerator delivers a reliable and powerful performance for every Dry Salt Therapy session. Its sleek lightweight design allows for easy installation for both small and large salt rooms. The HaloFX™ has numerous adjustable pre-set programs that make it customizable for any lung or skin condition. It is quiet, simple to operate and virtually maintenance free! Halotherapy Solutions Halogenerators are second to none. Maintenance and repair are minimal. Far and away the most user friendly, maintenance and repair free commercial / medical grade Halogenerator on the market.



## FEATURES

- Dimensions: 13"(W) x 21"(L) x 9"(D) | 14 lbs.
- Up to 440 sq. ft. room | Wall mount installation
- Low power consumption: 35-50W | 100-240V AC, 50/60 Hz
- Powdered feed rate: 0-300 g/h | Tank capacity: 200g
- 2 separate auxiliary controls for ventilation and lights
- Self-cleaning mill and non-rusting feeder
- Comprehensive training manual
- FREE Consultation on build-out

## LIST PRICE

- **\$6,500**

*\*For all customers, shipping cost is not included. For Non-US customers, you are responsible for customs/duties when applicable*

# HaloPod™

Introducing the HaloPod, one of our newer Halotherapy products, developed after much research with our wellness partners who requested a portable product that would easily fit into small spaces and that didn't require extensive or expensive room build-outs or modifications. Great for home use also!

## FEATURES

- Portable: Easy to move from room-to-room or to take outside of your spa to health fairs, senior centers to promote your business or add more revenue
- Includes portable pop-up pod to contain salt. The pod comes in 4 colors—black, red, navy blue, and royal blue
- Comfortable hanging attachment that maximizes dry salt absorption for individuals
- Durable Plastic vs. Stainless Steel (no clumping -easier to clean)

## LIST PRICE

- **\$4,900**



*\*For all customers, shipping cost is not included. For Non-US customers, you are responsible for customs/duties when applicable*